
Software

ACTIMED



Test- and management software

- $oxed{\square}$ management of devices, inspection management, costumer management
- oxdot scheduling of test activities
- ☑ DDE-link for controlling out of Microsoft application
- ☑ user specific printing with form generator
- $oxed{oxed}$ test examples included especially for Excel users

Test and measurement technic for medicine and industry







Technical Data

ACTIMED is a test- and management software for medical and industrial usage. It has been designed for in-house usage and mobile usage at your costumers (e.g. hospitals).

ACTIMED provides the creation, execution and mangagement of test procedures with devices from S.P.L. Electronik. You can control devices from other manufactures as long they are compatible with the communication gateway.



ACTIMED is based on GMEP and ACTIMED 1.0.



There are lots of costumer and device management features and a form generator to create specific prints. You can access to professional databases (e.g. InterBase) with the integrated ODBC-Database link.

(We reserve the right to make technical changes without prior notice 08/2011)