



Medical Pulp Macerator Consumables



EcoCleanse Macerator Chemical Options

Chemicals specifically designed for the unique operating requirements of your pulp or incontinence macerator, providing effective cleaning and deodorising.

EcoCleanse Ultimate chemicals are available in various product options, depending on the machine model and the specific market where it is installed. Please check with a DDC Dolphin representative for confirmation of the chemical product that can be supplied and what efficacy claims can be made.

Hygenex EcoCleanse Ultimate Chemicals



Product Description

EcoCleanse Ultimate chemical supplied in 5 litre bottles is available as a disinfectant, a deodoriser or as a cleaner.

EcoCleanse Ultimate Concentrate is supplied in a 1 litre cartridge compatible with the Pulpmatic Eco+ 'Click&Fit' system.

Order Code	Description	Box Qty	Pallet Qty
99-01-007	EcoCleanse Ultimate Disinfectant 5L	2	160
99-01-010	EcoCleanse Ultimate Deodoriser 5L	2	160
01-01-033	EcoCleanse Ultimate Cleaner 5L	2	160
99-01-080	EcoCleanse Ultimate Concentrate Disinfectant 1L	12	720

The Scalematic limescale inhibitor reduces the build-up of limescale throughout your Pulpmatic.

The Scalematic cartridge is fitted on the water supply prior to the machine connection point. By introducing polyphosphate into the water feed and sequestering the calcium in the water it helps prevent the formation of scale.

Replacing your Scalematic every 6 months will ensure your machinery stays in good working order for longer.

Scalematic



Product Description

The Scalematic limescale inhibitor safeguards your machinery against limescale build-up.

Order Code	Description	Qty
SCALEMATIC-SYS	Scalematic Limescale Inhibitor Unit	1
90-01-003	Scalematic Limescale Inhibitor replacement cartridge	1

Chemical Cabinets

Our wall mounted chemical cabinets are made of powder coated steel and are available with standard key lock or combination lock.

The cabinet is supplied with adaptor and flexible conduit to support the chemical supply tube linked directly to the machine.

Lockable Chemical Bottle Wall Cabinet



Product Description

Our wall mounted lockable chemical storage cabinets have the capacity for 2 x 5 litre bottles of EcoCleanse Ultimate or 12 x 1 litre bottles of EcoCleanse Ultimate Concentrate and are available with standard key lock or combination lock.

Order Code	Description	Qty
CAB-CHEM-001	Wall cabinet 5L lock & key	1
CAB-CHEM-002	Wall cabinet 5L Combination lock	1

Chemical Alarm Monitoring System (CAMS) takes the guess work out of when to replace the chemicals in a Pulpmatic Macerator, providing an audio and visual warning when the chemical bottle is empty.

The alarm can be set to either a loud or a softer pulsed beep, and there is also the option of a night mode which emits a red flashing light as a warning so as not to disturb sleep.

CAMS System



Product Description

Easy to fit the Chemical Alarm Monitoring System has the capacity to hold 2 x 5 litre bottles of EcoCleanse Ultimate and is available with a standard key lock or combination lock.

Order Code	Description	Qty
CAB-CHEM-005	CAMS Wall cabinet 5L lock & key	1
CAB-CHEM-006	CAMS Wall cabinet 5L Combination lock	1



Medical Pulp

Our medical pulp products are designed for general purpose use in surgeries, hospitals and care facilities.

We understand the importance of products that are not only cost effective, but which also give healthcare providers the confidence to perform tasks secure in the knowledge that safety and quality are key elements of the design and manufacturing processes allowing infection control strategies to be met.

Ergonomically designed for ease of use with a smooth, high-quality finish, helping to maintain hygiene standards and reduce the risk of cross contamination. Our products include items that are purpose built for managing a wide number of different, but entirely typical nursing situations that occur on an everyday basis, such as holding vomit, urine or faeces.

Our medical pulp products are made from recycled and maceratable pulp material in compliance with PAS29:1999 for healthcare disposable products and are Class 1 (According to Annex V111, Rule 1 of Medical Devices Regulation (EU) 2017/745). Consequently, the products specifically designed to contain human waste will do so for a minimum of 4 hours.

Bedpan (Figure 8) Liner & Accessories



Product Description

Suitable for general toileting of patients and can be used in bed or commode chair. The liner is placed inside the support before use. Ideal for patients confined to bed due to illness or after surgery.

Order Code	Description	Box Qty	Pallet Qty
91-01-005	Bed Pan 2 Litre Liner	100	20
91-02-008	Bedpan Support (figure 8)	1	N/A
91-02-003	Bedpan Dispenser Rack 2 Litre	1	N/A

Slipper Pan & Support (Midi)



Product Description

Recommended for general toileting purposes of patients with limited mobility in bed, can be slipped under the patient without rolling or lifting. Used along with a reusable plastic support these pulp liners improve hygiene and infection control.

Order Code	Description	Box Qty	Pallet Qty
91-01-003	Midi Slipper Pan 1.3 Litre Liner	100	17
91-02-009	Midi Slipper Pan Support	1	N/A

Slipper Pan Liner & Support (Maxi)



Product Description

Recommended for general toileting purposes of patients with limited mobility in bed, can be slipped under the patient without rolling or lifting. Used along with a reusable plastic support these pulp liners improve hygiene and infection control.

Order Code	Description	Box Qty	Pallet Qty
91-01-004	Maxi Slipper Pan 2 Litre Liner	100	16
91-01-010	Maxi Slipper Pan Support	1	N/A

Male Urinal Bottles & Accessories



Product Description

Enables independent patient toileting. The traditional round design is ergonomically shaped for use in bed, whilst the square design stands upright if required. Single urinal bottle holder also available, designed to stand upright or to be hooked on a bed and wall and floor mounted dispenser rack for medical pulp boxes.

Order Code	Description	Box Qty	Pallet Qty
91-01-001	Male Urinal Round 875 ml	100	5
91-01-002	Male Urinal Square 900 ml	120	6
WR-026	Single Urinal Bottle Holder	10	N/A
WR-025	Wall Mounted Frame	1	N/A

Commode Pan & Accessories



Product Description

Designed for general toileting in a commode chair. Also available Commode Pan Liner Dispenser rack.

Order Code	Description	Box Qty	Pallet Qty
91-01-013	Commode 1.7 Litre	200	10
91-02-005	Commode Pan Liner Dispenser Rack	1	N/A

General Purpose Bowls & Accessories



Product Description

General purpose bowls designed for ease of use available in 1 litre and 3 litre sizes. Also available General Purpose 1 litre dispenser rack.

Order Code	Description	Box Qty	Pallet Qty
91-01-006	General Purpose Bowl 1 Litre	200	24
91-01-012	General Purpose Bowl 3 Litre	100	20
91-02-006	General Purpose Dispenser Rack 1L	1	N/A



Detergent Proof Wash Bowl



Product Description

Holds warm water and detergent for up to 2 hours with a maximum 4 litre capacity.

Order Code	Description	Box Qty	Pallet Qty
91-01-011	Detergent Proof Wash bowl 4 Litre	100	12

Kidney Dish



Product Description

General use receptacle for holding dressings, wipes, surgical instruments and other medical items.

Order Code	Description	Box Qty	Pallet Qty
91-01-017	Kidney Dish 700ml	300	24

Graduated Jug



Product Description

General purpose jug with ergonomic handle for enhanced grip.

Order Code	Description	Box Qty	Pallet Qty
91-01-007	Jug 1 Litre	100	32

CS Receiver



Product Description

General use receptacle for holding dressings, wipes, surgical instruments and other medical items.

Order Code	Description	Box Qty	Pallet Qty
91-01-018	CS Receiver 900ml	300	12

Fluid Weighing Scales



Product Description

Hygenex fluid weighing scales. Electronic scales, quick and easy to use with easy clean surface. Maximum 5000 ml (for indication purposes only).

Order Code	Description	Qty
HYGSCALES	Fluid weighing scales	1

BioCommode Self Supporting Bedpan



Product Description

This single use, self supporting disposable bedpan is waterproof for up to 4 hours and supports up to 130kg without the need for plastic support.

Order Code	Description	Box Qty
90-02-011	BioCommode self supporting bedpan	70

Fluid Absorber Sachets

DDC Dolphin's Fluid Absorber sachets offer a convenient and effective way to handle harmful spills of liquid waste products.

Fluid absorber sachets have the ability to absorb and solidify spills into a gel, making disposal quick and easy. By using absorber sachets, the risk of handling infectious waste is reduced, while also eliminating unpleasant odours. The sachets dissolve upon contact with liquids and solidify them into a gel within a short period of time.

DDC Dolphin's liquid absorbers are made with a super absorbent polymer powder, allowing each sachet to absorb up to 100 times its own weight in liquid.

Fluid Absorber Sachets



Product Description

Fluid absorber sachets offer a high speed action to absorb and solidify harmful spills of liquid waste products into gel for ease of disposal.

Order Code	Description	Qty
91-01-108	Foil Bag containing 100 x 7 g sachets	1
91-01-109	Carton containing 10 foil bags of 100 x 7g sachets	1





**DDC Dolphin are your
infection control experts.**



DDC Dolphin Ltd

The Fulcrum,
Vantage Way,
Poole,
Dorset,
BH12 4NU



www.ddcdolphin.com



info@ddcdolphin.com



+44 (0) 1202 731555



Accreditations

